



aerospace
climate control
electromechanical
filtration
fluid & gas handling
hydraulics
pneumatics
process control
sealing & shielding



Filtration Housings

A guide to products and services



ENGINEERING YOUR SUCCESS.

Contents

General	
Precision Manufacturing	4-5
Project Management	6-7
A Dedicated Filter Range	8-9
Other Products from Parker domnick hunter	81
Parker Hannifin Group Product Range	82

Single Cartridge Housings

HSA	sanitary air / gas	12-13	HSI	in-line sanitary liquid	28-29
HSA⊕	sanitary air / gas	14-15	HSI⊕	in-line sanitary liquid	30-31
HBA	industrial air / gas	16-17	HIL	industrial liquid	32-33
HBA⊕	industrial air / gas	18-19	HIL⊕	industrial liquid	34-35
HSV	industrial vent	20-21	ZVP	industrial plastic	36
HSV⊕	industrial vent	22-23	Heating Jackets		37
HSL	sanitary liquid	24-25			
HSL⊕	sanitary liquid	26-27			

Small Scale Single Cartridge Housings					
HSA	sanitary air / gas	40-41	HSL	sanitary liquid	52-53
HSA⊕	sanitary air / gas	42-43	HSL⊕	sanitary liquid	54-55
HBA	industrial air / gas	44-45	HSI	in-line sanitary liquid	56-57
HBA⊕	industrial air / gas	46-47	HSI⊕	in-line sanitary liquid	58-59
HSV	industrial vent	48-49	HIF	industrial liquid	60-61
HSV⊕	industrial vent	50-51	HIF⊕	industrial liquid	62-63

Multi-Round Cartridge Housings					
ZVA	air / gas	66-67	VIL	industrial multi liquid	72-73
VIS	high flow steam	68-69	VSH	beverage	74-75
VSL	sanitary liquid	70-71			

Accessories	
Gauges	78
Valves	79
Spares	80
Certificates	80

Parker domnick hunter has a continuous policy of product development and although the Company reserves the right to change specification, it attempts to keep customers informed of any alterations. This publication is for general information only and customers are requested to contact our Process Filtration Sales Department for detailed information and advice on a product's suitability for specific applications. All products are sold subject to the company's Standard conditions of sale.

Precision Manufacturing

Experience and qualifications provide a flexible approach



Parker domnick hunter, Process Division manufacture stainless and carbon steel pressure vessels and filtration systems that are designed to meet International industry standards and specific customer application requirements.

A combination of highly skilled employees, dedicated manufacturing facility and 35 years experience of supplying process industries around the world Parker domnick hunter provide solutions that match your requirements for performance, quality and value.

Our fabrication facility manufacture a standard range of stainless steel housings to support our range of filters, which can be modified and adapted to meet any process requirements. Our strength is in providing a range of products that meet industry standards with a flexibility to meet your own process requirements.

Manufacturing Capability

- Pressure vessels from 0.1 to 10,000 litres
- Capacity: 5,000+ per year
- Automatic and hand welding techniques
- Assembly and hydro test facility
- Helium leak test, N.D.T., P.M.I. and stress relief
- Welding capability
 - manual / mechanical
 - MIG, MAG, TIG, MMA
 - micro plasma seam
 - keyhole plasma

Testing

- Helium leak test
- Surface finish
- Hydrostatic testing
- Pneumatic testing
- Ultrasonic testing
- Radiographic (x-ray)
- Swab testing
- Magnetic particle flow detection
- Riboflow testing
- Earth continuity testing

Manufacturing Best Practice


- ISO9001:2000
- ISO13485:2003
- ISO14001:2004

Vessels Built to Industry Standards

- PED (CE)
- EN / B445
- EN / 286
- BN / 1210
- ATEX
- ASME U
- ASME BPE

Stamp of Approval

- Certificate of Authorisation (U Stamp)
- National Board Certificate of Authorisation
- American Society of Mechanical Engineers



Project Management

Engineering your success



Parker domnick hunter, Process Division brings a wealth of experience in working on engineering projects around the world in partnership with some of the leading engineering, consultancy and project management groups. A highly trained workforce have the skills to match your exact requirements to the highest possible standards.

As part of the \$12 Billion turnover, Parker Hannifin Corporation, Parker domnick hunter can provide:

- Project management
- Process system design
- System fabrication
- Global support
- Operator training
- Dedicated technical support team
- Quality management systems

Our experience and expertise has seen us design and fabricate major systems for industries including:

- Pharmaceutical
- Chemical
- Food and beverage
- Industrial fermentation

A combination of hands on experience, design and manufacturing excellence have gained Parker domnick hunter a reputation for supplying high quality competitive filtration systems.

Leading Edge Design



Parker domnick hunter's, Sustaining Engineering Group are dedicated to providing a complete design service for coded pressure vessels, high integrity piping and all associated controls and instrumentation for project or contract work.

Using the latest in 3D CAD technology, Parker domnick hunter have the ability to support each project with:

- Visualisation - Photo Rendering
- Rapid Prototyping
- FEA - Finite Element Analysis
- CFD - Computational Flow Dynamics

Project Partnership

During the whole qualification phase of a new project, Parker domnick hunter provides formal validation plans, continuous support and assistance with all stages of qualification from factory acceptance through to site installation.



Dedicated Filter Range

Choice and flexibility to suit your application



Parker domnick hunter, Process Division manufacture a range of microfiltration cartridges for liquid and gas applications that utilise the latest production techniques, combining the most suitable membranes and filtration media with the latest easy to use formats.

All of Parker domnick hunter's filters meet strict validation guidelines providing a high degree of assurance that they will consistently achieve a high level of performance in a given application and meet the needs of the industry that they have been specifically designed for.

- Wide choice of filtration media and filter formats
- Technical and validation support
- Industry and application specific filters
- Fully retrofitable range of products
- Manufactured in state-of-the-art facilities

Scaleability provides flexibility
The ability to scale up from small area discs to final manufacturing with minimal revalidation is paramount.

Parker domnick hunter provides a wide range of filter formats to ensure that the transition from pilot scale through to full production is as smooth as possible.

Single use systems
Disposable systems can eliminate cleaning validation, reduce capital costs, minimise health & safety risks and lower the chance of product contamination. Single use systems also provide a more convenient way of processing a product.

Close working relationships
Parker domnick hunter have partnered engineering companies on large scale projects around the world that require filtration expertise and dedicated technical support.

Committed to process improvement
Direct access to our teams from new product development, laboratory services, technical support, manufacturing and quality provide the right solution delivered to you on time, every time.

Our goal is to continually improve your productivity, reduce your process costs and ensure the safety of your final product. Our Technical Support Group (TSG) made up from a multidisciplinary team of scientists and engineers working directly with your team to define your process needs and produce optimised solutions.



Single Cartridge Housings

5" to 40" cartridges



HSA - Sanitary air / gas housing

High specification air housing

HBA - Industrial air / gas housing

Specifically designed for the food & beverage industry

HSV - Vent housing

Flow efficient, self supportive sanitary housing

HSL - Sanitary liquid housing

Food, beverage & pharmaceutical finishes available

HSI - In-line sanitary liquid housing

Food, beverage & pharmaceutical finishes available

HIL - Industrial air / liquid housing

Ideal for water treatment & chemical applications



HSA Filter Housing

- sanitary air / gas

- Flow efficient sanitary range of air / gas housings
- Designed specifically for the food and beverage industry
- Sanitary vent, tri-clamp and drain connections as standard
- Sanitary tri-clamp body closure as standard



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	01K	011	012
Non Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)
PED Conformity Assessment Category			SEP	SEP	SEP
Volume [litres]			1.9	3.1	5.0

HSA Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl Head	Total	
5" (125 mm)	157	337	194	1.0	3.5	5.1
10" (250 mm)	157	487	313	1.6	3.5	5.7
20" (500 mm)	157	737	561	2.6	3.5	6.7

Dimensions are based on illustration shown (HSACE011YT-C-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSA			01				-		-		
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal	
CE	Standard	K	5" (125 mm)	Y	1 1/2" [38.1 mm]	T	Tri-Clamp	C	226	S	Silicone
		1	10" (250 mm)								
		2	20" (500 mm)								

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HSA® datasheet for more information.

HSA⊕ Filter Housing

- sanitary air / gas

- Sanitary range of air / gas housing
- Available in 4 different housing classes: Atex, CE, High Pressure and Oxygen Service
- Both beverage and pharmaceutical surface finishes available
- A choice of easy to use sanitary vent and drain options



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA, PTFE FDA, Silicone FDA, Viton FDA

Surface Finish

- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 µm Ra and Electropolished
 - External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration. All welds are crevice and undercut free. Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1. ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	135 °C (275 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	7.20 barg (104.42 psig)	5.60 barg (81.22 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I	CAT I	CAT I
Volume (litres)			1.9	3.1	5.0	6.9	8.8

CE Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	7.20 barg (104.42 psig)	5.60 barg (81.22 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I	CAT I	CAT I
Volume (litres)			1.9	3.1	5.0	6.9	8.8

High Pressure Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I	CAT I	CAT I
Volume (litres)			1.9	3.1	5.0	6.9	8.8



Oxygen Service Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	7.20 barg (104.42 psig)	5.60 barg (81.22 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I	CAT I	CAT I
Volume (litres)			1.9	3.1	5.0	6.9	8.8

HSA⊕ Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Typical Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
5" (125 mm)	229	337	205	0.9	1.9	3.3
10" (250 mm)	229	487	324	1.5	1.9	3.9
20" (500 mm)	229	737	572	2.5	1.9	4.9
30" (750 mm)	229	982	820	3.5	1.9	5.9
40" (1000 mm)	229	1232	1068	4.5	1.9	6.9

Dimensions shown are for a vessel with 2" tri-clamp ports, sample port and drain valve. For other formats, please contact Parker domnick hunter.



*Double bolted clamp required for HP and PTFE seal options

Ordering Information

HSA			01						-				-		-			
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal		Code Vent		Code Drain				
AT	ATEX	K	5" (125 mm)	Y	1 1/2" (38.1 mm)	T	Tri-Clamp	C	226	E	EPDM	C	Rectus 21 Vertical	H	Hosebarb			
CE	Standard	1	10" (250 mm)	C	2" (50.8 mm)					P*	PTFE	H	1 1/2" TCF & Hosebarb	R	Rectus 21			
HP*	High Pressure	2	20" (500 mm)							S	Silicone	I	1 1/2" TCF & Staubli RBE03	S	Staubli			
OX	Oxygen Service	3	30" (750 mm)							V	Viton	M	1 1/2" TCF & 1/2" TCF	R	RBE03			
		4	40" (1000 mm)									R	1 1/2" TCF & Rectus 21	T	1/2" TCF			
												T	Staubli RBE03 Vertical					
													1 1/2" TCF Only					
												Code Surface Finish		Internal	External			
												B		Beverage	0.4 µm	0.25 µm		
												P		Pharmaceutical	0.4 µm EP	0.25 µm		
																Code Tagged		
																T	Yes	
																X	No	

** Supplied complete with a double bolted clamp*

** Double bolted clamp required*

For Tagged Options customer identification numbers required at time of ordering

For Tagged Options customer identification numbers required at time of ordering

HBA Filter Housing

- industrial air / gas

- Flow efficient range of air / gas housings
- Designed to maximise flow and minimise pressure drop
- Designed specifically for the food and beverage industry



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA
- Vent / Drain Seals: PTFE

Surface Finish

- Internal: As Welded
 - External: Polished 0.8 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	01K	011	012
Non Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I
Volume (litres)			2.5	3.7	5.6

HBA Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
5" (125 mm)	231	396	223	0.9	1.9	3.6
10" (250 mm)	231	546	342	1.5	1.9	4.1
20" (500 mm)	231	796	590	2.5	1.9	5.2

Dimensions are based on illustration shown (HBACE011YB-C-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HBA			01				-		-		
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal	
CE	Standard	K	5" (125 mm)	Y	1 1/2" (38.1 mm)	B	BSPP	C	226	S	Silicone
		1	10" (250 mm)			N	NPT				
		2	20" (500 mm)								

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HBA® datasheet for more information.

HBA⊕ Filter Housing

- industrial air / gas

- Flow efficient range of air / gas housings
- Available in 4 different housing classes: Atex, CE, High Pressure and Oxygen Service
- Beverage, pharmaceutical and industrial surface finishes available
- A number of inlet / outlet port connections
- Wide range of vent and drain options



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA, PTFE FDA, Silicone FDA, Viton FDA

Surface Finish

- Industrial Finish
 - Internal: As Welded
 - External: Pickled & Passivated
- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 µm Ra and Electropolished
 - External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration. All welds are crevice and undercut free. Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1. ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	135 °C (275 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	6.60 barg (95.72 psig)	5.30 barg (76.87 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I	CAT I
Volume (litres)			2.5	3.7	5.6	7.5	9.4

CE Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	6.60 barg (95.72 psig)	5.30 barg (76.87 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I	CAT I
Volume (litres)			2.5	3.7	5.6	7.5	9.4

High Pressure Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I	CAT I
Volume (litres)			2.5	3.7	5.6	7.5	9.4



Oxygen Service Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Dangerous	Gas / Vapour	150 °C (302 °F)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	8.00 barg (116.03 psig)	6.60 barg (95.72 psig)	5.30 barg (76.87 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I	CAT I
Volume (litres)			2.5	3.7	5.6	7.5	9.4

HBA⊕ Filter Housings

Physical Characteristics

Bowl Height		Dimensions (mm)			Typical Weight (Kg)		
		'A'	'B'	'C'	Bowl	Head	Total
5"	(125 mm)	259	398	223	1.0	5.4	7.0
10"	(250 mm)	259	548	342	1.6	5.4	7.6
20"	(500 mm)	259	798	590	2.6	5.4	8.6
30"	(750 mm)	259	1043	838	3.6	5.4	9.6
40"	(1000 mm)	259	1293	1068	4.6	5.4	10.6

Dimensions shown are for a vessel with 1 1/2" BS4504 DIN2633 ports, 1/4" BSP vent and drain. For other formats, please contact Parker domnick hunter.



*Double bolted clamp required for HP and PTFE seal options

Ordering Information

HBA											01						-			-		-	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal		Code Vent			Code Drain								
AT	ATEX	K	5" (125 mm)	Y	1 1/2" (38.1 mm)	B	BSPP (F)	C	226	E	EPDM	B	1/4" BSPP (F)	N	G 1/2" BSPP								
CE	Standard	1	10" (250 mm)	C	2" (50.8 mm)	D	DIN11851(M)			P	PTFE	C	Rectus 21 Vertical ⁽²⁾	N	G 1/2" NPT								
HP*	High Pressure	2	20" (500 mm)			F	ANSI RF 150 ⁽¹⁾			S	Silicone	H	1 1/2" TCF & Hosebarb ⁽²⁾	H	Hosebarb								
OX	Oxygen Service	3	30" (750 mm)			H	ANSI RF 300			V	Viton	I	1 1/2" TCF & Staubli RBE03 ⁽²⁾	R	Rectus 21								
										* Double bolted clamp required		M	1 1/2" TCF & 1/2" TCF ⁽²⁾	S	Staubli RBE03								
												N	1/2" NPT (F)	R	Staubli RBE03								
										(2) Not available in Industrial Finish.		N	1 1/2" TCF & Rectus 21 ⁽²⁾	T	1/2" TCF								
												S	Staubli RBE03 Vertical										
												T	1 1/2" TCF Only ⁽²⁾										
												Code Surface Finish		Internal		External							
												B	Beverage	0.4 µm		0.25 µm							
												I	Industrial	As Welded		0.8 µm							
												P	Pharmaceutical	0.4 µm EP		0.25 µm							
												Code Tagged											
												T	Yes										
												X	No										
												For Tagged Options customer identification numbers required at time of ordering											

HSV Filter Housing

- vent applications

- Vent housings
- Direct connection to tank boss allows housing to be self-supportive
- Corrosion resistant 316L stainless steel
- Easy assembly and maintenance



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.8 µm Ra
 - External: As welded
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration. All welds are crevice and undercut free. Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Recommended Operation Guidelines Sizing

Sizing vent vessels particularly for vacuum sensitive tanks can require specialist advice. It is important that VENT housings are sized on maximum gas flow capacity under actual operation conditions.

Vacuum Protection



Where a tank is vacuum sensitive, there is a risk of tank collapse. In such cases the fitting of an appropriately rated bursting disc (or similar) and, if necessary a pressure relief valve, is highly recommended.

HSV Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
5" (125 mm)	132	242	194	0.9	0.7	2.2
10" (250 mm)	132	392	313	1.5	0.7	2.8
20" (500 mm)	132	642	561	2.5	0.7	3.8

Dimensions are based on illustration shown (HSVDH011YT-C-S). For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSV		01				-		-	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge	
DH	Vent Housing	K	5" (125 mm)	Y	1 1/2" (38.1 mm)	T	Tri-Clamp	C	226
		1	10" (250 mm)					S	Silicone
		2	20" (500 mm)						

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range. Please see HSV® datasheet for more information.

HSV⊕ Filter Housing

- industrial vent



- Industrial vent housings
- Available in Atex version
- Beverage, pharmaceutical and industrial surface finishes available
- Available in various connection types

Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Note: Seal used only to position bowl clamp arrangement.

Surface Finish

- Industrial Finish
 - Internal: As Welded
Pickled & Passivated
 - External: Polished 0.8 µm Ra
- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 µm Ra and Electropolished
 - External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Basis

ATEX 94/9/EC (where applicable)



Working Condition PED 97/23/EC			Volume (litres)				
Variant	State	Temperature	01K	011	012	013	014
Standard	Gas / Vapour	150 °C [302 °F]	1.6	2.8	4.7	6.6	8.5
ATEX	Gas / Vapour	135 °C [275 °F]	1.6	2.8	4.7	6.6	8.5

HSV⊕ Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Typical Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
5" (125 mm)	132	274	194	0.9	2.6	4.0
10" (250 mm)	132	424	313	1.5	2.6	4.6
20" (500 mm)	132	674	561	2.5	2.6	5.6
30" (750 mm)	132	919	809	3.5	2.6	6.6
40" (1000 mm)	132	1169	1057	4.5	2.6	7.6

Dimensions shown are for a vessel with 1½" BS4504 DIN 2633 outlet port.
For other formats, please contact Parker domnick hunter.



Ordering Information

HSV		01						-		-							
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal		Code Surface Finish		Internal	External	Code Tagged	
DH	Standard	K	5" (125 mm)	C	2" (50.8 mm)	B	BSPP (F)	C	226	E	EPDM	B	Beverage	0.4 µm	0.25 µm	T	Yes
AT	ATEX	1	10" (250 mm)			D	DIN11851(M)			P	PTFE	I	Industrial	As Welded	0.8 µm	X	No
		2	20" (500 mm)	Y	1½" (38.1 mm)	F	ANSI RF150			S	Silicone	P	Pharmaceutical	0.4 µm EP	0.25 µm		
		3	30" (750 mm)			L	BS4504			V	Viton						
		4	40" (1000 mm)				DIN2633										
						N	NPT (F)										
						T	Tri-Clamp										
						W	BS / ISO Pipe										

For Tagged Options customer identification numbers required at time of ordering

For Tagged Options customer identification numbers required at time of ordering

HSL Filter Housing

- sanitary liquid

- Single element sanitary liquid housing
- Designed specifically for the food and beverage and pharmaceutical industry
- Sanitary vent and tri-clamp connections as standard
- Sanitary tri-clamp body closure as standard



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	011	012	013
Non Dangerous	Liquid / Gas	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid / Gas	150 °C (302 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I
Volume [litres]			2.9	4.8	6.7

HSL Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
10" (250 mm)	156	417	313	1.0	1.5	3.8
20" (500 mm)	156	667	561	1.6	1.5	4.8
30" (750 mm)	156	912	809	2.6	1.5	5.7

Dimensions are based on illustration shown (HSLCE011BT-C-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSL		01				-		-	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge	
CE	Standard	1	10" (250 mm)	B	1" (25.4 mm)	T	Tri-Clamp	C	226
		2	20" (500 mm)					S	Silicone
		3	30" (750 mm)						

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HSL® datasheet for more information.

HSI Filter Housing

- in-line sanitary liquid

- In-line sanitary liquid housing
- High quality crevice free construction
- Sanitary body closure as standard



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	011	012	013
Non Dangerous	Liquid / Gas	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid / Gas	150 °C (302 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I
Volume [litres]			2.8	4.7	6.6

HSI Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl Head	Total	
10" (250 mm)	132	449	313	1.0	1.5	3.6
20" (500 mm)	132	699	561	1.6	1.5	4.6
30" (750 mm)	132	944	809	2.6	1.5	5.6

Dimensions are based on illustration shown (HSICE011YT-C-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSI		01				-		-	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge	
CE	Standard	1	10" (250 mm)	Y	1 1/2" (38.1 mm)	T	Tri-Clamp	C	226
		2	20" (500 mm)					S	Silicone
		3	30" (750 mm)						

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HSI® datasheet for more information.

HSI⊕ Filter Housing

- in-line sanitary liquid

- In-line sanitary liquid housing
- Available in 3 different housing classes: ATEX, CE and High Pressure
- Both beverage and pharmaceutical surface finishes available
- Sampling and drain port options



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Surface Finish

- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra

- Pharmaceutical Finish

- Internal: Polished 0.4 µm Ra and Electropolished
- External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration. All welds are crevice and undercut free. Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.
ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	135 °C (275 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
Non Dangerous	Liquid	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	135 °C (275 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	SEP	SEP	CAT I	CAT I
Volume (litres)			1.6	2.8	4.7	6.6	8.5



CE Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C (302 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
Non Dangerous	Liquid	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	150 °C (302 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	SEP	SEP	CAT I	CAT I
Volume (litres)			1.6	2.8	4.7	6.6	8.5

High Pressure Working Condition PED 97/23/EC			Maximum Pressure				
Fluid Group	State	Temperature	01K	011	012	013	014
Non Dangerous	Gas / Vapour Liquid	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category			SEP	SEP	CAT I	CAT I	CAT I
Volume (litres)			1.6	2.8	4.7	6.6	8.5

Physical Characteristics

Bowl Height	Dimensions (mm)	Typical Weight (Kg)		
		'A'	'B'	'C'
5" (125 mm)	132 327 194	0.9	1.5	2.9
10" (250 mm)	132 477 313	1.5	1.5	3.5
20" (500 mm)	132 727 561	2.5	1.5	4.5
30" (750 mm)	132 972 809	3.5	1.5	5.5
40" (1000 mm)	132 1222 1057	4.5	1.5	6.5

Dimensions shown are for a vessel with a vent sample port. For other formats, please contact Parker domnick hunter.



*Double bolted clamp required for HP and PTFE seal options

Ordering Information

HSI 01 -									
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge	
AT	ATEX	K	5" (125 mm)	Y	1 1/2" (38.1 mm)	T	Tri-Clamp	C	226
CE	Standard	1	10" (250 mm)						
HP*	High Pressure	2	20" (500 mm)						
		3	30" (750 mm)						
		4	40" (1000 mm)						
Code Seal		Code Vent		Code Drain		Code Surface Finish		Internal	
E	EPDM	H	1 1/2" TCF & Hosebarb	H	Hosebarb	B	Beverage	0.4 µm	0.25 µm
P*	PTFE	I	1 1/2" TCF & Staubli RBE03	R	Rectus 21	P	Pharmaceutical	0.4 µm EP	0.25 µm
S	Silicone	M	1 1/2" TCF & 1/2" TCF	S	Staubli				
V	Viton	R	1 1/2" TCF & Rectus 21	T	RBE03				
		X	No Vent	X	1/2" TCF				
					No Drain				
Code Tagged		Code Tagged		Code Tagged		Code Tagged		Code Tagged	
T	Yes	T	Yes	T	Yes	T	Yes	T	Yes
X	No	X	No	X	No	X	No	X	No

For Tagged Options customer identification numbers required at time of ordering

HIL Filter Housing

- industrial liquid

- Industrial single element liquid housing
- 1" BSPP or NPT inlet / outlet standard connections
- Suitable replacement for plastic housings
- Suitable for cartridge types DOE or 222



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA

Surface Finish

- Internal: As Welded
 - External: Polished 0.8 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	011	012	013
Non Dangerous	Liquid / Gas	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid / Gas	150 °C (302 °F)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I
Volume [litres]			3.2	5.1	7.0

HIL Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Weight (Kg)		
	'A'	'B'	'C'	Bowl	Head	Total
10" (250 mm)	150	441	297	1.5	1.2	3.8
20" (500 mm)	150	691	550	2.5	1.2	4.9
30" (750 mm)	150	936	814	3.5	1.2	6.0

Dimensions are based on illustration shown (HILCE011BB-B-E).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HIL	01					
Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seal	
CE Standard	1 10" (250 mm) 2 20" (500 mm) 3 30" (750 mm)	B 1" (25.4 mm)	B N BSPP NPT	B D DOE 222	E EPDM	

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HIL datasheet for more information.

HIL⊕ Filter Housing

- industrial liquid



- Industrial single element liquid housing
- Available in 3 different housing classes:
Atex, CE and High Pressure
- Industrial and industrial-electropolished
surface finishes available
- Suitable for cartridge types DOE or 222
- Fabricated 'C' style version available
(Not Cast Head)

Specification

Materials of Construction

- Housing: 316L Stainless Steel (Cast Head)
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Surface Finish

Two Finished Available:

- **Industrial Finish**
Head-Cast, Pickled & Passivated
Bowl Internal: As Welded
Bowl External: Pickled & Passivated
- **Industrial Electropolished Finish**
Head-Cast, Pickled & Passivated
Bowl Internal: Electropolished
Bowl External: Polished 0.8 µm Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.
ATEX 94/9/EC (where applicable)



Working Condition PED 97/23/EC			Maximum Pressure			
Fluid Group	State	Temperature	011	012	013	014
Non Dangerous	Gas / Vapour	135 °C [275 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	135 °C [275 °F]	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
Non Dangerous	Liquid	135 °C [275 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	135 °C [275 °F]	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I
Volume (litres)			3.2	5.1	7.0	8.9

Working Condition PED 97/23/EC			Maximum Pressure			
Fluid Group	State	Temperature	011	012	013	014
Non Dangerous	Gas / Vapour	150 °C [302 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Gas / Vapour	150 °C [302 °F]	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
Non Dangerous	Liquid	150 °C [302 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	150 °C [302 °F]	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)	5.00 barg (72.51 psig)
PED Conformity Assessment Category			SEP	CAT I	CAT I	CAT I
Volume (litres)			3.2	5.1	7.0	8.9

High Pressure Working Condition PED 97/23/EC			Maximum Pressure			
Fluid Group	State	Temperature	011	012	013	014
Non Dangerous	Gas / Vapour	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
Non Dangerous	Liquid	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category			CAT I	CAT I	CAT I	CAT I
Volume (litres)			3.2	5.1	7.0	8.9

	Bowl Height	Dimensions [mm]			Typical Weight (Kg)		
		'A'	'B'	'C'	Bowl	Head	Total
10"	[250 mm]	230	441	297	1.5	5.4	7.8
20"	[500 mm]	230	691	550	2.5	5.4	8.9
30"	[750 mm]	230	936	814	3.5	5.4	10.0
40"	[1000 mm]	230	1186	1058	4.5	5.4	11.1

Dimensions shown are for a vessel with 1" BS4504 DIN 2633 ports,
1/4" BSP vent and drain. For other formats, please contact
Parker domnick hunter.



*Double bolted clamp required for HP and PTFE seal options

Ordering Information

HIL

01

-

-

-

Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal		Code Vent		Code Drain	
AT	ATEX	1	10" (250 mm)	B	1" (25.4 mm)	B*	BSPP (F)	B	DOE	E	EPDM	B	1/4" BSPP	B	1/4" BSPP
CE	Standard	2	20" (500 mm)			F	ANSI RF150 ⁽¹⁾	C	226 (Fabricated Head)	P*	PTFE	N	1/4" NPT	N	1/4" NPT
HP*	High Pressure	3	30" (750 mm)	Y	1 1/2" (38.1 mm)	H	ANSI RF 300	D	222	S	Silicone	X*	No Vent		
		4	40" (1000 mm)			L	BS4504 DIN2633			V	Viton				
						N* NPT (F) T*o Tri-clamp						*Only available with 'C' Style cartridge configuration			
										* Double bolted clamp required					

ZVP Housings

- industrial plastic

- Single cartridge polypropylene / nylon housing
- Accepts DOE filters with knife edge sealing
- Accepts plug-in cartridges with positive o-ring seals
- Meets water conditioning foundation standards for hydraulic leak test and ultimate burst pressure
- Cost-effective filtration of liquids for pharmaceutical, chemical and beverage applications




Available Options

Type	Material	Design Pressure (barg)	Design Temperature
ZVP-1	Reinforced Polypropylene	8.6	51.7 °C [125 °F]
ZVP-2	Reinforced Polypropylene	8.6	51.7 °C [125 °F]
ZVP-3	Polycarbonate / Reinforced Polypropylene	8.6	51.7 °C [125 °F]
ZVP-4	Nylon	8.5	71.1 °C [160 °F]
ZVP-5	Reinforced Polypropylene	8.6	51.7 °F [125 °F]
ZVP-7	Reinforced Polypropylene	8.6	51.7 °C [125 °F]
ZVP-10	Reinforced Polypropylene	8.6	51.7 °C [125 °F]
ZVP-11	Pure Polypropylene	8.9	38.0 °C [100.4 °F]

Type	Crossport (A) (mm)	Overall Height (B) (mm)	Approx Weight (kg)	Approx Volume (L)	Connection Size / Type	Vent Button
ZVP-1	112.7	180	0.6	0.7	3/8" BSPP	Yes
ZVP-2	130.0	311	1.5	1.6	1/2" BSPP	Yes
ZVP-3	130.0	321	1.2	1.6	3/4" BSPP	Yes
ZVP-4	130.0	305	1.8	1.6	1/2" BSPP	No
ZVP-5	130.0	569	1.9	2.6	3/4" BSPP	Yes
ZVP-7	130.0	311	1.5	1.6	3/4" BSPP	Yes
ZVP-10	130.0	569	1.9	2.6	3/4" BSPP	Yes
ZVP-11	130.0	316	1.4	1.6	3/4" BSPP	No*


* 1/2" BSPP vent and drain sockets with o-rings and plugs



Heating Jackets

The design of heating systems for vent filters working in high humidity environments is often treated as an afterthought, but the correct operation of these filters can be critical to many processes. Applications such as the venting of Water For Injection (WFI) holding tanks rely on a well designed heated housing to prevent condensation build up, filter blockage and the risk of microbial growth. It is also one of the key design requirements highlighted in current FDA recommendations. Heating may also be required during hot water sanitisation and CIP to prevent excess differential pressure being generated from high levels of bulk condensate.

- Heating system for vent applications
- Waterproof protection to IP65
- Fully insulated 'cool touch' outer surface
- Accurate temperature control using PT100RT6



Specification - Heating Jacket

Materials of Construction

- Jacket Material
 - Silicone: Silicone Rubber
 - Glass Silk: PTFE Coated Glass Silk
- Insulation Material
 - Silicone: Silicone Foam
 - Glass Silk: PTFE Coated Glass Silk

Maximum Withstand Temperature De-Energised
200 °C [392 °F]

Temperature Sensor
PT100

Thermal Cut-Out Temperature Setting
150 °C ± 5 °C [302 °F ± 41 °F]

Test Voltage
1500 V

Insulation Value
Greater than 100 mΩ

Protection Rating

Silicone: IP65
Glass Silk: Not Applicable
Inter-Connection Plugs: IP67

Design Standards

EN 60519-1 and EN 60519-2

Operating Voltage
110 V or 230 V

Power Output
5" (125 mm): 63 W
10" (250 mm): 279 W
20" (500 mm): 558 W
30" (750 mm): 837 W

Specification - Temperature Control Unit

Materials of Construction

- Material: Polycarbonate

Operating Voltage
110 V or 230 V

Maximum Withstand Temperature of Controller
55 °C [131 °F]

Maximum Continuous Current Out
7 A

Over Current Protection @ 230 °C Ambient
4 seconds @ 12 A, 1 second @ 24 A

Set Temperature Display
8 mm Red LED Display

Actual Temperature Display
10 mm Green LED Display

PID (Proportional Integral Derivative) Control
Via autotune parameters (set by user)

Protection Rating
IP65

Design Standards
EN 61010-1

Alarm Contacts (Normally Open) When Applicable

Switching Capacity Max.
250 VAC 0.5 A (load resistance)
125 VAC 1 A (load resistance)
60 VAC 1 A (load resistance)



Small Scale Single Housings

Demi 2.5" & 5" small scale filter housings



HSA - Sanitary air / gas housing

Pharmaceutical & beverage grade finishes

HBA - Industrial air / gas housing

Specifically designed for the food & beverage industry

HSV - Vent housing

Pharmaceutical & beverage grade finishes

HSL - Sanitary liquid housing

Pharmaceutical & beverage grade finishes

HSI - In-line sanitary liquid housing

Pharmaceutical & beverage grade finishes

HIF - Industrial air / liquid housing

Industrial grade finish as standard



Demi HSA Filter Housing

- sanitary air / gas

- Flow efficient sanitary range of air / gas housings
- Designed specifically for the food and beverage industry
- Sanitary tri-clamp, vent and drain connections as standard
- Sanitary tri-clamp body closure as standard



Specification

Materials of Construction
■ Housing: 316L Stainless Steel
■ Seals: Silicone FDA

Surface Finish
■ Internal: Polished 0.4 µm Ra
■ External: Polished 0.25 µm Ra
All finishes pickled & passivated.

Welding
All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification
Supplied as standard with vessel inspection certificate.

Material Test Certification
EN10204 3.1 supplied upon request.

Design Code
Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis
ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume [litres]			0.75	0.50

Demi HSA Filter Housings

Physical Characteristics

	Bowl Height	Dimensions (mm)			Typical Weight (Kg)
		'A'	'B'	'C'	
A Size	5" (125 mm)	152	227	130	1.3
B Size	2 1/2" (65 mm)	152	172	70	1.2

Dimensions are based on illustration shown (HSACE01ABT-T-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSA		<div></div>	01	<div></div>	<div></div>	<div></div>	-	<div></div>	-	<div></div>	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge		Code Seal	
CE	Standard	A	5" (125 mm)	B	1"	T	Tri-Clamp	T	126	S	Silicone
		B	2 1/2" (65 mm)								

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range. Please see HSA® datasheet for more information.

Demi HSA⊕ Filter Housing

- sanitary air / gas



- Sanitary range of air / gas housing
- Available in 4 different housing classes: Atex, CE, High Pressure and Oxygen Service
- Both beverage and pharmaceutical surface finishes available
- A choice of easy to use sanitary vent and drain options

Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Surface Finish Options

- Beverage Finish
 - Internal: Polished 0.4 μ m Ra
 - External: Polished 0.25 μ m Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 μ m Ra and Electropolished
 - External: Polished 0.25 μ m Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.
ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50



CE Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	150 °C [302 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.5

High Pressure Working Condition PED 97/23/EC				Maximum Pressure	
Fluid Group		State	Temperature	01A	01B
Non Dangerous		Gas / Vapour	205 °C [401 °F]	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category				SEP	SEP
Volume (litres)				0.75	0.5


Oxygen Service Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.5

Bowl Height	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
A Size 5" (125 mm)	152	340	130	1.3
B Size 2 1/2" (65 mm)	152	285	70	1.2

Dimensions are based on illustration shown (HSACE01ABTTE-HH-P-X)
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information



Demi HBA Filter Housing

- industrial and beverage air / gas

- Flow efficient range of air / gas housings
- Designed to maximise flow and minimise pressure drop
- Designed specifically for the food and beverage industry



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Body Seal: Silicone FDA
- Vent / Drain Seal: PTFE FDA

Surface Finish

- Internal: Unpolished 1 µm Ra Typical
 - External: Polished 0.8 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	150 °C [302 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume [litres]			0.75	0.50

Demi HBA Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
A Size 5" (125 mm)	175	300	130	1.5
B Size 2 1/2" (65 mm)	175	245	70	1.4

Dimensions are based on illustration shown (HBACE01ABB-TS).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HBA		01				-		-	
Code Vessel Class		Code Length (Nominal)		Code Connection Size		Code Standard		Code Cartridge	
CE	Standard	A	5" (125 mm)	B	1"	B	BSPP	T	126
		B	2 1/2" (65 mm)			N	NPT	S	Silicone

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HBA® datasheet for more information.

Demi HBA⊕ Filter Housing

- industrial and beverage air / gas

- Flow efficient range of air / gas housings
- Available in 4 different housing classes: Atex, CE, High Pressure and Oxygen Service
- Beverage, pharmaceutical and industrial surface finishes available
- A number of inlet / outlet port connections
- Wide range of vent and drain options



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Body Seal: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA
- Vent / Drain Seal: PTFE FDA

Surface Finish Options

- Industrial Finish
 - Internal: As Welded
Pickled & Passivated
 - External: Polished 0.8 µm Ra
- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 µm Ra
and Electropolished
 - External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.
ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

CE Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

High Pressure Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous	Gas / Vapour	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50



Oxygen Service Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Dangerous	Gas / Vapour	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

Demi HBA⊕ Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)	'A'	'B'	'C'	Typical Weight (Kg)
A Size	5" (125 mm)	220	300	130	1.5
B Size	2 1/2" (65 mm)	220	245	70	1.4

Dimensions are based on illustration shown (HBACE01ABFTE-BB-P-X).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HBA 01 -							
Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seal	Code Vent	Code Drain
AT ATEX CE Standard HP High Pressure OX Oxygen Service	A 5" (125 mm) B 2 1/2" (65 mm)	B 1" T 3/4" A 1/2" X 3/8" Q 1/4"	B BSPP (F) D ⁽²⁾ DIN11851 F ⁽³⁾⁽⁴⁾ ANSI RF 150 H ⁽⁵⁾ ANSI RF 300 L ⁽³⁾ BS4504 DIN2633 M* SMS Pipe [3008] N NPT (F) T ⁽³⁾ Tri-Clamp W ISO / BS Pipe	T 126	E EPDM P* PTFE S Silicone V Viton * Double bolted clamp required	B 1/4" BSPP (F) C Rectus 21 Vertical H ⁽¹⁾ 1 1/2" TCF & Hosebarb I ⁽¹⁾ 1 1/2" TCF & Staubli RBE03 M ⁽¹⁾ 1 1/2" TCF & 1/2" TCF N 1/4" NPT (F) R ⁽¹⁾ 1 1/2" TCF & Rectus 21 S Staubli RBE03 Vertical T 1 1/2" TCF Only	B 1/4" BSPP N 1/4" NPT H Hosebarb R Rectus 21 S Staubli RBE03 T 1/2" TCF
Code Surface Finish		Internal		External		Code Tagged	
B Beverage		0.4 µm		0.25 µm		T Yes X No	
I Industrial		As Welded		0.8 µm			
P Pharmaceutical		0.4 µm EP		0.25 µm			
<div>For Tagged Options customer identification numbers required at time of ordering</div>							

For Tagged Options customer identification numbers required at time of ordering

Demi HSV Filter Housing

• vent housing

- Direct connection to tank boss allows housing to be self-supportive
- Corrosion resistant 316L stainless steel
- Easy assembly and maintenance



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.8 µm Ra
 - External: As welded
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration. All welds are crevice and undercut free. Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Working Condition PED 97/23/EC		Volume (litres)	
State	Temperature	01A	01B
Gas / Vapour	150 °C [302 °F]	0.75	0.50

Recommended Operation Guidelines Sizing

Sizing vent vessels particularly for vacuum sensitive tanks can require specialist advice. It is important that VENT housings are sized on maximum gas flow capacity under actual operation conditions.

Vacuum Protection



Where a tank is vacuum sensitive, there is a risk of tank collapse. In such cases the fitting of an appropriately rated bursting disc (or similar) and, if necessary a pressure relief valve, is highly recommended.

Demi HSV Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)		Typical Weight (Kg)
	'B'	'C'	
A Size 5" [125 mm]	195	130	1.0
B Size 2 1/2" [65 mm]	140	70	0.9

Dimensions are based on illustration shown [HSV01ABT-T-S]. For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSV	01			-		-	
Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seat		
DH Vent Housing	A 5" [125 mm] B 2 1/2" [65 mm]	B 1"	T Tri-Clamp	T 126	S Silicone		

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range. Please see HSV® datasheet for more information.

Demi HSV⊕ Filter Housing

- vent housing

- Available in ATEX version
- Beverage, pharmaceutical and industrial surface finishes available
- Available in various connection types



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Note: Seal used only to position bowl clamp arrangement.

Surface Finish Options

- Industrial Finish
 - Internal: As Welded
Pickled & Passivated
 - External: Polished 0.8 µm Ra
- Beverage Finish
 - Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 µm Ra and Electropolished
 - External: Polished 0.25 µm Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Basis

ATEX 94/9/EC (where applicable)



Working Condition PED 97/23/EC			Volume (litres)	
Variant	State	Temperature	01A	01B
Standard	Gas / Vapour	150 °C [302 °F]	0.75	0.5
ATEX	Gas / Vapour	135 °C [275 °F]	0.75	0.5

Demi HSV⊕ Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)		Typical Weight (Kg)
	'B'	'C'	
A Size 5" [125 mm]	203	130	1.7
B Size 2 1/2" [65 mm]	148	70	1.6

Dimensions are based on illustration shown [HSV01ABFTE-P-X].
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSV			01							-		-			
Code Vessel Class		Code Length [Nominal]		Code Connection Size		Code Standard		Code Cartridge		Code Seal		Code Surface Finish		Internal	External
DH	Standard	A	5" [125 mm]	B	1"	B	BSPP (F)	T	126	E	EPDM	B	Beverage	0.4 µm	0.25 µm
AT	ATEX	B	2 1/2" [65 mm]			D	DIN11851			P	PTFE	I	Industrial	As Welded	0.8 µm
						F	ANSI RF150			S	Silicone	P	Pharmaceutical	0.4 µm EP	0.25 µm
						L	BS4504			V	Viton				
						N	DIN2633								
						T	NPT (F)								
							Tri-Clamp								
														Code Tagged	
														T	Yes
														X	No

For Tagged Options customer identification numbers required at time of ordering

Demi HSL Filter Housing

- sanitary liquid

- Single element sanitary liquid housing
- Sanitary tri-clamp, vent and drain connections as standard
- Sanitary tri-clamp body closure as standard



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone FDA

Surface Finish

- Internal: Polished 0.4 µm Ra
 - External: Polished 0.25 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour & Liquid	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

Demi HSL Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
A Size 5" (125 mm)	117	191	130	2.0
B Size 2 1/2" (65 mm)	117	136	70	1.8

Dimensions are based on illustration shown (HSLCE01ABT-T-S).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HSL		01				-		-	
Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seat				
CE Standard	A 5" (125 mm) B 2 1/2" (65 mm)	B 1"	T Tri-Clamp	T 216	S Silicone				

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HSL® datasheet for more information.

Demi HSL⊕ Filter Housing

- sanitary liquid



- Single element sanitary liquid housings
- Available in 3 different housing classes: Atex, CE and High Pressure
- Both beverage and pharmaceutical surface finishes available
- Wide range of connection, vent and drain options available

Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA

Surface Finish Options

- Beverage Finish
 - Internal: Polished 0.4 μ m Ra
 - External: Polished 0.25 μ m Ra
- Pharmaceutical Finish
 - Internal: Polished 0.4 μ m Ra and Electropolished
 - External: Polished 0.25 μ m Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.
ATEX 94/9/EC (where applicable)



ATEX Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour & Liquid	135 °C [275 °F]	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

CE Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour & Liquid	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

High Pressure Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous	Gas / Vapour Liquid	205 °C (401 °F)	16.00 barg [232.06 psig]	16.00 barg [232.06 psig]
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

	Bowl Height	Dimensions (mm)			Typical Weight (Kg)
		'A'	'B'	'C'	
A Size	5" (125 mm)	132	290	130	2.3
B Size	2 1/2" (65 mm)	132	235	70	2.2

*Dimensions are based on illustration shown (HSLCE01ABDTE-HH-P-X).
For accurate dimensions, please contact Parker domnick hunter.*



Ordering Information

[illegible]

Demi HSI Filter Housing

- in-line sanitary liquid

- In-line sanitary liquid housing
- High quality crevice free construction
- Sanitary body closure as standard



Specification

Materials of Construction
■ Housing: 316L Stainless Steel
■ Seals: Silicone FDA

Surface Finish
■ Internal: Polished 0.4 µm Ra
■ External: Polished 0.25 µm Ra
All finishes pickled & passivated.

Welding
All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification
Supplied as standard with vessel inspection certificate.

Material Test Certification
EN10204 3.1 supplied upon request.

Design Code
Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis
ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Liquid / Gas	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

Demi HSI Filter Housings

Physical Characteristics

	Bowl Height	Dimensions (mm)		Typical Weight (Kg)
		'B'	'C'	
A Size	5" (125 mm)	207	130	1.0
B Size	2½" (65 mm)	152	70	0.8

*Dimensions are based on illustration shown (HSICED1AYT-T-S).
For accurate dimensions, please contact Parker domnick hunter.*



Ordering Information

HSI		<div></div>	01	<div></div>	<div></div>		<div></div>	-	<div></div>	-	<div></div>	
Code Vessel Class		Code Length (Nominal)			Code Connection Size		Code Standard		Code Cartridge		Code Seal	
CE	Standard	A	5" (125 mm)	B	2½" (65 mm)	Y	1½"	T	Tri-Clamp	T	216	S Silicone

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HSI® datasheet for more information.

Demi HIF Filter Housing

- industrial air / liquid

- Industrial single element air / liquid housing
- 1/2" BSPP or NPT inlet / outlet standard connections
- Suitable replacement for plastic housings
- Suitable for Parker domnick hunter 'Z' style 116 'O' rings



Specification

Materials of Construction

- Housing: 316L Stainless Steel (Cast Head)
- Body Seal: EPDM FDA
- Vent / Drain Seal: PTFE FDA

Surface Finish

- Internal: Unpolished 1 µm Typical
 - External: Polished 0.8 µm Ra
- All finishes pickled & passivated.

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Certification

Supplied as standard with vessel inspection certificate.

Material Test Certification

EN10204 3.1 supplied upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK Statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.



Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Liquid / Gas	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

Demi HIF Filter Housings

Physical Characteristics

Bowl Height	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
A Size 5" (125 mm)	75	226	130	1.5
B Size 2 1/2" (65 mm)	75	171	70	1.4

Dimensions are based on illustration shown (HIFCE01AAB-Z-E).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

HIF	01			-		-	
Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seat		
CE Standard	A 5" (125 mm) B 2 1/2" (65 mm)	A 1/2"	B BSPP N NPT	Z 116	E EPDM		

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

For additional features, Parker domnick hunter offer this housing as part of its Standard PLUS Range.
Please see HIF® datasheet for more information.

Demi HIF⊕ Filter Housing

- industrial air / liquid



- Industrial single element air / liquid housing
- Available in 3 different housing classes:
Atex, CE and High Pressure
- Industrial and industrial-electropolished
surface finishes available
- Suitable for Parker domnick hunter
'Z' style 116 'O' rings

Specification

Materials of Construction

- Housing: 316L Stainless Steel
(Cast Head)
- Body Seal: EPDM FDA
PTFE FDA
Silicone FDA
Viton FDA
- Vent / Drain Seal: PTFE FDA

Surface Finish Options
Two Finishes Available:

- **Industrial Finish**
 Head-Cast, Pickled & Passivated
 Bowl Internal: As Welded
 Pickled & Passivated
 Bowl External: Polished 0.8 μ m Ra
- **Industrial Electropolished Finish**
 Head-Cast, Pickled, Passivated
 & Electropolished
 Bowl Internal: Electropolished
 Bowl External: Polished 0.8 μ m Ra

Welding

All assembly welds are full penetration.
All welds are crevice and undercut free.
Weld finish & detail drawings available upon request.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis



ASME VIII Division 1.
ATEX 94/9/EC (where applicable)

ATEX Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour & Liquid	135 °C (275 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50
CE Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous & Dangerous	Gas / Vapour & Liquid	150 °C (302 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50
High Pressure Working Condition PED 97/23/EC			Maximum Pressure	
Fluid Group	State	Temperature	01A	01B
Non Dangerous	Gas / Vapour	205 °C (401 °F)	16.00 barg (232.06 psig)	16.00 barg (232.06psig)
PED Conformity Assessment Category			SEP	SEP
Volume (litres)			0.75	0.50

Physical Characteristics

Bowl Height			Dimensions (mm)			Typical Weight (Kg)
			'A'	'B'	'C'	
A Size	5"	(125 mm)	75	226	130	1.5
B Size	2 1/2"	(65 mm)	75	171	70	1.4

Dimensions are based on illustration shown (HIFCE01AABZS-BB-E-X).
For accurate dimensions, please contact Parker domnick hunter.



Ordering Information

Diagram illustrating the 10-PORT MANIFOLD configuration with various options. The manifold is shown with 10 ports, each corresponding to a specific code or option. The options are defined in the tables below.

Code Vessel Class	Code Length (Nominal)	Code Connection Size	Code Standard	Code Cartridge	Code Seal	Code Vent	Code Drain
AT ATEX	A 5" (125 mm)	A 1/2"	B BSPP (F)	Z 116	E EPDM	B 1/4" BSPP	B 1/4" BSPP
CE Standard	B 2 1/2" (65 mm)	X 3/8"	N NPT (F)		P* PTFE	N 1/4" NPT	N 1/4" NPT
HP High Pressure					S Silicone		
					V Viton		

* Double bolted clamp required

Code Surface Finish	Internal	External
E Industrial Electropolished	Electropolished	Linished 0.8 µm
I Industrial	As Welded	Linished 0.8 µm

Code Tagged
T Yes
X No

For Tagged Options customer identification numbers required at time of ordering

Multi Housings

3 to 30 multi round cartridge housings



ZVA - Sanitary range air / gas housing

Specifically designed for the pharmaceutical industry

VSL - Multi-element sanitary liquid housing

Designed specifically for the pharmaceutical industry

VIL - Multi-element industrial liquid housing

General purpose industrial housing

VSH - Multi-element liquid housing

Designed for prefiltration & clarification applications

VIS - High flow steam

Specifically designed for steam filtration



ZVA Housings

- air / gas

- Sanitary range air / gas housing
- Specifically designed for the Pharmaceutical industry
- Laboratory and pilot scale to large industrial applications
- Flow efficient design with low pressure drop
- Steam jacketed and electrically heated options



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM

Surface Finish

- Multis - Basic Specification
 - Internal: Inside of outlet assy and distribution box to be mechanically polished 0.8 µm Ra. Immerse entire vessel to achieve 100% pickle and passivation.
 - External: Grit blast 5 µm Ra mean

- Multis - Full Specification
 - Internal: Electropolish 0.6 µm Ra
 - External: Bright Polished 0.4 µm Ra

Maximum Allowable Working Pressure (MAWP) PS
6 barg (87.0 psig)

Maximum Allowable Working Temperature (MAWT) TS
120 °C [248 °F]

Maximum Allowable Working Pressure Steam
3 barg (43.5 psig) @ 144 °C (291 °F)

Total Volume (litres)

031	032	033
31.0	40.0	49.0
051	052	053
45.0	58.0	72.0

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001. PED / PER conformity assessments based on Fluid Group 2 Gas (harmless) including steam. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII Division 1.



Custom Design

Parker domnick hunter offers a specialist and fabrication service allowing individual customer system specifications to be met.

Physical Characteristics

Bowl Height	Dimensions (mm)	Typical Weight (Kg)
10" (250 mm)	336 794 320	43.0
20" (500 mm)	336 1044 570	47.0
30" (750 mm)	336 1294 820	50.0

Dimensions shown are for a ZVACE 3 round, 3" ANSI inlet / outlet connections. For the full range of dimensions and weights, please contact Parker domnick hunter.



Ordering Information

Code	Vessel Class	Code	N° of Cartridges	Code	Length (Nominal)	Code	Connection Size	Code	Connection Type	Code	Connection Standard	Code	Vent / Drain Conn. Type	Code	Surface Finish
CE	Standard	03	3	1	10" (250 mm)	D	3"	F	Flanged	A	NPT / ANSI			C	Economy Spec
		05	5	2	20" (500 mm)	E	4"	W	Weld Prepared	N	Nominal Bore	B	BSPP	F	Full Spec
		09	9	3	30" (750 mm)	G	6"			I	ISO				
		15	15			H	8"								

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

- Specifically designed to maximise flow rates and minimise pressure drop
- Compatible with JUMBO element to maximise steam capacity

Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM

Surface Finish

- Internal: Inside of outlet and distribution box to be mechanically mirror polished 0.8 µm Ra. Immerse vessel to achieve 100% pickle and passivation.
- External: Grit blast 5 µm Ra mean

Maximum Allowable Working Pressure (MAWP) PS
7 barg (101.5 psig)


Maximum Allowable Working Temperature (MAWT) TS
170.5 °C (339 °F)

Design Code
Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory pressure equipment regulations (PER) 1999 N° 2001. PED / PER
Conformity assessments based on Fluid Group 2 Gas (harmless) allowing for in-situ steam sterilisation. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis
ASME VIII Division 1.

VIS Housings



- high flow steam



Physical Characteristics

Type	Dimensions (mm)			Typical Weight [Kg]
	'A'	'B'	'C'	
VISCE-01J-DD	300	763	310	30.0
VISCE-01J-ED	330	895	140	50.0
VISCE-3J-DD	515	1049	410	100.0
VISCE-3J-ED	700	1237	490	150.0

For the full range of dimensions and weights, please contact Parker domnick hunter.



Ordering Information

VIS			-			-												
Code Vessel Class		Code N° of Cartridges		Code Length (Nominal)		Code Connection Size		Code Connection Type										
CE	Standard	01	1	J	Jumbo	D*	3"	D	BS4504 PN16 Flange									
		03	3			E*	4"			A	ANSI cl. 150 Flange							
						G**	6"											
						H**	8"											
* Single housings only																		
** Round housings only																		

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

VSL Housings

- sanitary liquid



- Multi-element sanitary liquid housing
- Designed specifically for the pharmaceutical industry
- Electropolished internal finish

Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM

Surface Finish

- Internal: Electropolished 0.4 μm Ra
- External: Polished 0.25 μm Ra

Economy Spec

An economy version is available with a lower specification, external lished to 0.8 μm Ra.

Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	3 / 5 Round	3 / 5 Round	3 / 5 Round
Non Dangerous	Liquids	80 °C (176 °F)	7.00 barg (101.50 psig)	7.00 barg (101.50 psig)	7.00 barg (101.50 psig)
Dangerous / Non Dangerous	Gas / Vapour	144 °C (297 °F)	3.00 barg (43.50 psig)	3.00 barg (43.50 psig)	3.00 barg (43.50 psig)
PED Conformity Assessment Category			CAT I	CAT II	CAT II
Volume (litres)			10.7	18.5	26.3

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory pressure equipment regulations (PER) 1999 N° 2001. PED / PER Conformity assessments based on Fluid Group 2 Gas (harmless) allowing for in-situ steam sterilisation. Only housings over PS.V 50 bar / litres bear the CE mark.

Design Basis

ASME VIII.

Custom Design



Parker domnick hunter offers a specialist and fabrication service allowing individual customer system specifications to be met.

VSL Housings

Physical Characteristics

Type	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
10" (250 mm)	271	493	320	18.0
20" (500 mm)	271	743	570	22.0
30" (750 mm)	271	993	820	26.0

Dimensions shown are for a 3 Round VSL, 2" TCF inlet / outlet connections. For the full range of dimensions and weights, please contact Parker domnick hunter.



Ordering Information

VSL				-				-							
Code Vessel Class		Code N° of Cartridges		Code Length (Nominal)		Code Connection Size		Code Connection Type		Code Connection Standard		Code Connection Type		Code Seal	
CE Standard		03 3 05 5		1 10" (250 mm) 2 20" (500 mm) 3 30" (750 mm)		C 2"		T Tri-Clamp		B British Standard D DIN 07* Economy Series		B Sanitary Bleed T Tri-Clamp		E EPDM P PTFE S Silicone V Viton	
										*Code finishes here, no mention of drain and seals					

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

VIL Multi Filter Housing

- industrial multi liquid

- Multi-element industrial liquid housing
- Laboratory and pilot scale to large industrial applications
- Flow efficient design with low pressure drop



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: EPDM

Surface Finish

- Internal / External:
 - DOE Economy: As fabricated then pickled to remove weld discolouration
 - DOE Standard: As fabricated then electropolished
 - P-7 (226) o-ring: As fabricated then electropolished

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory Pressure Equipment Regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.

Custom Design

Parker domnick hunter offers a specialist and fabrication service allowing individual customer system specifications to be met.



Working Condition PED 97/23/EC			Maximum Pressure		
Fluid Group	State	Temperature	3 / 5 Round	8 Round	12 Round
Non Dangerous	Liquids	80 °C (176 °F)	10.00 barg (145.00 psig)	8.50 barg (123.00 psig)	6.00 barg (87.00 psig)
PED Conformity Assessment Category			SEP	SEP	SEP
Volume (litres)			031	032	033
			14.0	19.0	24.0
			051	052	053
			20.0	28.0	36.0
			082	083	084
			70.0	87.0	104.0
			122	123	124
			98.0	122.0	146.0

VIL Multi Filter Housings

Physical Characteristics

Type		Dimensions (mm)			Typical Weight [Kg]
		'A'	'B'	'C'	
10"	(250 mm)	222	720	70	14.0
20"	(500 mm)	222	970	70	16.0
30"	(750 mm)	222	1220	70	18.0
40"	(1000 mm)	222	1470	70	20.0

Dimensions shown are for a 3 Round VIL, 2" BSPP inlet / outlet connections. For the full range of dimensions and weights, please contact Parker domnick hunter.



Ordering Information

VIL									
Code Vessel Class		Code N° of Cartridges		Code Length [Nominal]		Code Connection Size		Code Connection Type	
CE	Standard	03	3	1	10" (250 mm)	C	2"	B	BSPP Female
		05	5	2	20" (500 mm)	D	3"	F	Flanged
		08	8	3	30" (750 mm)			N	NPT Female
		12	12	4	40" (1000 mm)	E	4"		
Code Connection Style		Code Seals		Code Surface Finish		Code Flange Type*		Code Cartridge Length [dh standard]	
B	DOE	E	EPDM	E	Electropolished	15	ANSI cl. 150	S	Slimline L
C	P-7	P***	PTFE	P**	Pickled	45	BS4504		10", 20", 30", 40"
		S	Silicone						
		V	Viton						


Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

- Note:
- * Only applicable for housings with flange connection type
 - ** Pickled finish option only available for DOE housings
 - *** PTFE seal option requires flanged closure

VSH Multi Housings

- beverage

- Multi-element sanitary liquid housing
- Designed specifically for the food and beverage industry
- High quality crevice free construction
- Available in 3 to 30 round versions
- Steam sterilisable



Specification

Materials of Construction

- Housing: 316L Stainless Steel
- Seals: Silicone

Surface Finish

- Internal: Mechanically Polished Ra <0.8 µm
- External: Mechanically Polished

Steam Sterilisation

Refer to Parker domnick hunter for individual housing parameters.

Design Code

Housings designed in accordance with the European Council Pressure Equipment Directive (PED) 97/23/EC and the UK statutory pressure equipment regulations (PER) 1999 N° 2001.

Design Basis

ASME VIII Division 1.

Custom Design

Parker domnick hunter offers a specialist and fabrication service allowing individual customer system specifications to be met.

Note: For 12, 18, 24 and 30 Round options, please contact Parker domnick hunter for detailed technical drawings.



Working Condition PED 97/23/EC			Maximum Pressure			
Fluid Group	State	Temperature	031	032	033	034
Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Dangerous	Gas / Vapour	0 - 150 °C (0 - 302 °F)	4.80 barg (69.62 psig)	3.90 barg (56.56 psig)	2.80 barg (40.61 psig)	2.10 barg (30.45 psig)
Non Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Non Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Non Dangerous	Gas / Vapour	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Non Dangerous	Gas / Vapour	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Volume (litres)			7.3	12.6	17.8	23.1
Fluid Group	State	Temperature	051	052	053	054
Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Dangerous	Gas / Vapour	0 - 150 °C (0 - 302 °F)	4.50 barg (65.26 psig)	2.40 barg (34.80 psig)	1.70 barg (24.65 psig)	1.30 barg (18.85 psig)
Non Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Non Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Non Dangerous	Gas / Vapour	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	9.90 barg (143.58 psig)	6.80 barg (98.62 psig)	5.20 barg (75.41 psig)
Non Dangerous	Gas / Vapour	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	6.80 barg (98.62 psig)	5.20 barg (75.41 psig)
Volume (litres)			11.0	20.0	29.1	38.2
Fluid Group	State	Temperature	081	082	083	084
Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Dangerous	Gas / Vapour	0 - 150 °C (0 - 302 °F)	2.30 barg (33.35 psig)	1.40 barg (20.30 psig)	1.00 barg (14.50 psig)	0.70 barg (10.15 psig)
Non Dangerous	Liquid	0 - 40 °C (0 - 104 °F)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)	10.00 barg (145.03 psig)
Non Dangerous	Liquid	150 °C (302 °F)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)	7.50 barg (72.51 psig)
Non Dangerous	Gas / Vapour	0 - 40 °C (0 - 104 °F)	9.40 barg (136.33 psig)	5.60 barg (81.22 psig)	4.00 barg (58.01 psig)	3.10 barg (44.96 psig)
Non Dangerous	Gas / Vapour	150 °C (302 °F)	7.50 barg (72.51 psig)	5.60 barg (81.22 psig)	4.00 barg (58.01 psig)	3.10 barg (44.96 psig)
Volume (litres)			21.3	35.3	49.7	63.9
PED Conformity Assessment Category			CAT I	CAT I	CAT I	CAT I

VSH Multi Housings

Physical Characteristics

Type	Dimensions (mm)			Typical Weight (Kg)
	'A'	'B'	'C'	
10" (250 mm)	606	840	290	27.0
20" (500 mm)	606	1060	540	30.0
30" (750 mm)	606	1310	790	33.0
40" (1000 mm)	606	1560	1040	36.0

Dimensions shown are for an 8 Round VSH, DN40 DIN2633 inlet / outlet connections. For the full range of dimensions and weights, please contact Parker domnick hunter.



Ordering Information

VSH - - - - -


Code Vessel Class	Code N° of Cartridges	Code Length (Nominal)	Code Connection Size	Code Connection Type	Code Connection Standard	Code Vent / Drain Conn. Type	Code Seal
CE Standard	03 3 05 5 08 8 12 12 18 18 24 24 30 30	1 10" (250 mm) 2 20" (500 mm) 3 30" (750 mm) 4 40" (1000 mm)	B 1" C 1.5" & 2" X 2.5" D 3" E 4"	D DIN Male F Flanged M SMS Male R RJT Male T Tri-Clamp W Weld Prepared	A NPT / ANSI B British D DIN I ISO	T Tri-Clamp	E EPDM S Silicone V Viton

Note: 3, 5 and 8 Round housings available with 10", 20", 30" and 40" bowls
12, 18, 24 and 30 Round housings available with 20" and 30" bowls only

N° of Cartridges Connection Size Availability	1.5"	2"	2.5"	3"	4"
03	✓				
05	✓	✓			
08	✓	✓	✓		
12				✓	✓
18				✓	✓
24				✓	✓
30				✓	✓

Note: For accessories, i.e. gauges, please contact Parker domnick hunter - Process Division for full availability.

Accessories



Gauges

Industry & product specific options

Valves

Comprehensive range of manual valves

Spares

Replacement parts & accessories

Certificates

Comprehensive range of supporting documentation



Parker domnick hunter provide a comprehensive range of pressure gauges to support their standard air / gas and liquid housings.

HIL Pressure Gauge			
Type	Connection	Pressure	Ordering Code
All stainless steel wetted parts with glycerine fill fluid. Design temperature 135 °C (275 °F). Cooling tower required for temperatures up to 205 °C (401 °F). See spares page 80.	1/4" BSP	0 - 10 barg	XPGSS03BS17
		0 - 16 barg	XPGSS03BS18
	1/4" NPT	0 - 10 barg	XPGSS03NP2
		0 - 16 barg	XPGSS03NP3

HBA & HIF Pressure Gauge			
Type	Connection	Pressure	Ordering Code
All stainless steel wetted parts with glycerine fill fluid (includes adapter for plain 1/2" BSPP connection). Design temperature 135 °C (275 °F). Cooling tower required for temperatures up to 205 °C (401 °F). See spares page 80.	1/2" BSP	0 - 10 barg	XPGSS03BS15
		0 - 16 barg	XPGSS03BS16
	1/2" NPT	0 - 10 barg	XPGSS03NP2
		0 - 16 barg	XPGSS03NP3

HSL Single & VSH Multi Pressure Gauge		
Type	Pressure / Temperature	Ordering Code
Sanitary gauge with double sanitary valves and sight glass for beverage applications. Glycerine gauge fill fluid.	0 - 10 barg / 150 °C (302 °F)	XGSSS08TC1

HSL Single & VSH Multi Pressure Gauge Valve Manifold		
Type	Pressure / Temperature	Ordering Code
Sanitary gauge manifold with double sanitary valves without sight glass, for applications where glass not allowed. (Manifold and quick release only. Does not include gauge, clamp and seal.)	0 - 10 barg / 150 °C (302 °F)	XMISS06TC1

HSL Single & VSH Multi Pressure Gauge		
Type	Pressure / Temperature	Ordering Code
Sanitary gauge with single sanitary valve. Glycerine gauge fill fluid.	0 - 10 barg / 150 °C (302 °F)	XPGSS08TC1

Single, VSH & VSL Multi Pressure Gauge		
Type	Pressure / Temperature	Ordering Code
Sanitary gauge with 1" and 1 1/2" tri-clamp connection and hygienic diaphragm. All stainless steel wetted parts with KN92 FDA approved fill fluid.	0 - 10 barg / 150 °C (302 °F)	XPGSS06TC4



Parker domnick hunter provide a comprehensive range of manual valves to support their standard air / gas and liquid housings.

Industrial 1 Piece Ball Valve		
Type	Connection	Ordering Code
316 stainless steel 1 piece ball valve with PTFE ball. Male / female.	1/4" BSP	XVASS03BS
	1/4" NPT	XVASS03NP

2 Piece Ball Valve		
Type	Connection	Ordering Code
316 stainless steel semi-sanitary ball valve. (for use on up-stream connection on sanitary liquid beverage or pharmaceutical housings)	1 / 1 1/2" Tri-Clamp - 20 mm Hosebarb	XVASS06TC

Butterfly Valves		
Type	Variant	Ordering Code
Stainless steel butterfly valve with silicone seals and polymer handle.	1" OD x 1.6 Weld End	XVASS050D1/VHPL
	1 1/2" OD x 1.6 Weld End	XVASS060D1/VHPL
	2" OD x 1.6 Weld End	XVASS070D1/VHPL
	2 1/2" OD x 1.6 Weld End	XVASS080D1/VHPL
	3" OD x 1.6 Weld End	XVASS090D1/VHPL
	DN25 DIN11851	XVASS050N1/VHPL
	DN40 DIN11851	XVASS060N1/VHPL
	DN50 DIN11851	XVASS070N1/VHPL
	DN65 DIN11851	XVASS080N1/VHPL
	DN80 DIN11851	XVASS090N1/VHPL
	1" Tri-clamp Ferrule	XVASS05TC4/VHPL
	1 1/2" Tri-clamp Ferrule	XVASS06TC4/VHPL
	2" Tri-clamp Ferrule	XVASS07TC4/VHPL
	2 1/2" Tri-clamp Ferrule	XVASS08TC4/VHPL
	3" Tri-clamp Ferrule	XVASS09TC4/VHPL

Sanitary Bleed Valve			
Type	Seals	Variant	Ordering Code
316 stainless steel sanitary bleed valve with Neoprene grip. Available with EPDM, Silicone, Viton or Perlast Seals. Available with Rectus 21, Staubli RBE03 or 8 mm hosebarb.	EPDM	Staubli RBE03 Male	XVASS30NA1
		Rectus 21 Male	XVASS30RT
		8 mm Hosebarb	XVASS30HB
	Silicone	Staubli RBE03 Male	XVASS30ST1
		Rectus 21 Male	XVASS30RT1
		8 mm Hosebarb	XVASS30HB1
	Viton	Staubli RBE03 Male	XVASS30NA4
		Rectus 21 Male	XVASS30RT2
		8 mm Hosebarb	XVASS30HB2
	Perlast	Staubli RBE03 Male	XVASS30NA2
		Rectus 21 Male	XVASS30RT3
		8 mm Hosebarb	XVASS28SL15

Sample Valve		
Type	Connection	Ordering Code
316 stainless steel sanitary valve with 1" / 1 1/2" tri-clamp connection and 12 mm hosebarb. For use on down-stream connection on sanitary liquid housings.	1 / 1 1/2" Tri-Clamp - Stepped 12 mm Hosebarb	XVASS05TC3

Gemu Diaphragm Valve			
Type	Connection	Variant	Ordering Code
316 stainless steel sanitary diaphragm valve with 1/2" (miniclap) tri-clamp connection and silicone or EPDM diaphragm.	1/2" (miniclap) tri-clamp	Silicone	XVASS04TC6
		EPDM	XVASS04TC1
		Viton	XVASS04TC7
		PTFE	XVASS04TC8



Parker domnick hunter provide a comprehensive range of spare parts to support their standard air / gas and liquid housings.

4" Spares		2 1/2" & 4" Spares		3 Round VSH Spares	
Size & Type	Part Code	Size & Type	Part Code	Size & Type	Part Code
4" Single Pin Tri-Clamp 4" Double Bolt Tri-Clamp	XTCSS10SL XTCSS10HP15	1/4" BSP Plug 1/4" NPT Plug	XPLSS03BS4 XPLSS03NP1	3 Round Body 'V' Clamp	XBCSS51BL
4" TCF Gasket EPDM 4" TCF Gasket Silicone 4" TCF Gasket Viton 4" Gasket PTFE	XTSEP10SA XGKSI3004 XTSV10SL XTSPT10SL	1/4" BSP PTFE Plug Seal ATEX Earth Kit (Replacement)	XGKPT03BP XEKSS00AT	3 Round Spider Plate	XSPSS51BL
HIL 222 Spring HIL DOE Nut	XSNS0070D XNTSS011L	1 1/2" Tri-Clamp Blanking Kit - EPDM 1 1/2" Tri-Clamp Blanking Kit - Silicone 1 1/2" Tri-Clamp Blanking Kit - Viton 1 1/2" Tri-Clamp Blanking Kit - PTFE	XAKSS06TC6 XAKSS06TC3 XAKSS06TC7 XAKSS06TC8	Body O-Ring BS362 - Silicone 1 1/2" Head to Elbow Gasket - PTFE	XORSI12BL XGKPT06BL
2 1/2" Spares		1/2" Tri-Clamp Blanking Kit - EPDM 1/2" Tri-Clamp Blanking Kit - Silicone 1/2" Tri-Clamp Blanking Kit - Viton 1/2" Tri-Clamp Blanking Kit - PTFE	XAKSS04TC1 XAKSS04TC2 XAKSS04TC3 XAKSS04TC4	5 Round VSH Spares	
Size & Type	Part Code	Size & Type	Part Code	Size & Type	Part Code
2 1/2" Single Pin Tri-Clamp 2 1/2" Double Bolt Tri-Clamp	XTCSS08SA XTCSS08HP	1 1/2" Tri-Clamp Gasket - EPDM 1 1/2" Tri-Clamp Gasket - Silicone 1 1/2" Tri-Clamp Gasket - Viton 1 1/2" Tri-Clamp Gasket - PTFE	XGKEPTC XTSSI06 XGKVI06TC2 XGKPT06TC	5 Round Body 'V' Clamp	XBCSS52BL
2 1/2" TCF Gasket EPDM 2 1/2" TCF Gasket Silicone 2 1/2" TCF Gasket Viton 2 1/2" Gasket PTFE	XGKEP08NA XGKSI08 XGKVI08SA XGKPT10SA	1 1/2" Tri-Clamp Gasket - EPDM 1 1/2" Tri-Clamp Gasket - Silicone 1 1/2" Tri-Clamp Gasket - Viton 1 1/2" Tri-Clamp Gasket - PTFE	XGKEP04TC1 XGKSI04TC XGKVI04TC XGKPT04TC	5 Round Spider Plate	XSCSS52BL
		Sanitary Seal Kit - EPDM Sanitary Seal Kit - Silicone Sanitary Seal Kit - Viton Sanitary Seal Kit - Perlast	XOREP30 XORSI30 XORVI30 XORPE30NA1	Body O-Ring BS370 - Silicone 2" Head to Elbow Gasket - PTFE	XORSI12BL1 XGKPT07
		1 1/2" TC Blank 1 1/2" TC Clamp	XTBSS05TC XTCSS05TC	8 Round VSH Spares	
		1 1/2" TC Blank 1 1/2" TC Clamp	XTBSS04TC XTCSS04TC	Size & Type	Part Code
		Cooling Tower 1/4" BSPP Cooling Tower 1/4" NPT	XCTSS03BS XCTSS03NP	8 Round Body 'V' Clamp	XBCSS52BL
		VSH Spares		8 Round Spider Plate	XSPSS52BL
		Size & Type	Part Code	Body O-Ring BS378 - Silicone 3" Head to Elbow Gasket - PTFE	XORSI13BL XGKPT07
		3, 5 & 8 Round Vent and Drain Clamp / Gasket Kit	XAKSS06TC5		

Parker domnick hunter provide a comprehensive range of certificates to support their standard air / gas and liquid housings.

Certificate	
Type	Ordering Code
Vessel Inspection Certificate (included with vessel)	66 950 0013
Vessel Inspection Certificate (replacement)	66 950 0013
Material Certification Pack (EN10204 3.1)	66 950 0014
Certificate of Conformity	66 950 0015
Passivation Report	66 950 0016
Cleanliness Certificate	66 950 0017
Surface Finish Certificate	66 950 0018
Weld Procedure Certificate Pack	66 950 0019
Quality Plan	66 950 0026
Replacement IOMI (Installation, Operation and Maintenance Instructions)	17 950 0769



Industrial Products

Parker domnick hunter, Industrial Division, is a well established global business capable of meeting the compressed air treatment product needs of all industries. Our commitment to customer satisfaction goes beyond initial supply and installation. Comprehensive after sale support includes servicing, spare parts, quality testing and technical advice.

Bespoke design services are also available for customised projects to ensure customer specifications are met. Services are delivered locally by our global network of qualified service engineers.



LAB GAS GENERATORS
Hydrogen, nitrogen & zero air
The range of analytical gas generators from Parker domnick hunter includes UHP hydrogen, nitrogen and zero air and enables users to produce a cost-effective, continuous supply of premium quality gas from a compact, on-site source.

- Increases safety with the elimination of high pressure gas storage or cylinder handling
- Cost-effective due to low life-cycle ownership
- UHP hydrogen generators facilitate optimised analysis
- Convenient, on-demand gas supply



MAXIGAS
Nitrogen gas generators
Produces on-site nitrogen gas from compressed air and is the cost-effective alternative to traditional nitrogen sources for multiple applications. Excellent energy efficiency and a low life-cycle ownership cost facilitate considerable cost savings of up to 90%.

- Low life-cycle ownership cost and elimination of costs associated with a cylinder supply
- On-demand functionality limits waste
- Energy efficient; operates from a small compressor



MIXED GAS DISPENSERS
CO₂ & nitrogen
Designed to provide bar owners with the ideal supply of mixed gas blends of CO₂ and nitrogen for beer dispensing. The system uses a nitrogen generator which, when connected to CO₂ cylinders, can produce mixed blends of CO₂ and nitrogen in a number of predetermined ratios.

- Improved quality and economy
- Nitrogen purity of 99.8%
- A more efficient operation
- Improved shelf life



WS WATER SEPARATORS
Bulk liquid removal
Providing efficient bulk liquid removal at all flow conditions, OIL-X EVOLUTION WS Water Separators also minimise energy consumption and help reduce your carbon footprint.

- Tested in accordance with ISO8573.9
- Performance independently verified
- Low pressure loss / low operational cost



OIL-X EVOLUTION
Compressed air filters
Providing air quality that meets or exceeds the requirements of ISO8573-1, the international standard for compressed air quality, OIL-X EVOLUTION is also the most energy efficient compressed air filter in the world, helping to reduce your carbon footprint.

- The most energy efficient filters available
- High quality ISO8573.1:2001 compressed air
- Running costs that start low and stay low



CRD
Refrigeration dryers
Avoid corrosion, machinery failure and product spoilage by removing water from any compressed system at affordable prices. The CRD range provides the very latest in drying technology and is suitable for all compressor types.

- Clean, dry compressed air, stops damage and corrosion
- Environmentally friendly R407C refrigerant
- Energy efficient, low running costs



PNEUDRI
Desiccant dryers
Providing water vapour removal in accordance with Classes 1, 2 & 3 of ISO8573-1 the international standard for compressed air quality, PNEUDRI modular compressed air dryers offer unrivalled performance, flexibility and expandability in a unique space saving design. Low operational costs and integrated energy management systems also ensure energy consumption is kept to a minimum.

- Highest quality air
- Totally stops corrosion / damage
- Low installation costs
- Energy efficient



ES2000 SERIES
Oil / water separators
Providing a legal and responsible way to dispose of oil contaminated compressor condensate, ES2000 series oil water separators are a cost effective alternative to expensive waste disposable companies.

- Help to protect and maintain the environment
- Efficiently separate oil and water on-site and return up to 99.9% of the condensate to foul sewers
- Meet trade effluent discharge regulations
- Rapid payback over conventional disposal methods



BREATHING AIR PURIFIERS
Breathable air
Providing breathable quality compressed air in compliance with international standards, breathing air purifiers supply effective protection from harmful substances, maintaining employee health.

- High efficiency coalescing filter, for removal of oil / water
- Adsorption bed of activated carbon, for removal of oil vapour and odours
- Catalytic element, for removal of carbon monoxide



NBC FILTRATION
Biological & chemical protection
The need to protect key personnel from attacks by chemical and biological weapons has never been greater. Given the escalation of this type of threat from terrorist groups and unstable nations, the development of the NBC filtration system provides effective protection.

- Fully regenerative
- Increased capacity
- Compact modular design



HYPERCHILL
Precision chilled water
Hyperchill maximises productivity and minimises costs, as well as easy conformity to regulations on water quality. Hyperchill is the perfect solution to industrial chilled water needs.

- Increases productivity, reduces costs
- Adaptable to individual customer needs



PCO₂
Carbon dioxide polishing filter
Providing quality incident protection for beverage grade carbon dioxide, PCO₂ offers protection against carbon dioxide contamination and impurities of up to 10 times the allowable levels.

- Ensures compliance with quality guidelines published by the International Society for Beverage Technologies (ISBT)
- Protects drinks manufacturing processes from vapour impurities

For further information on the full range of Industrial and Gas Generation products available, please contact Parker domnick hunter Industrial Division

tel: +44 (0)191 402 9000
fax: +44 (0)191 482 6296
email: dhindsales@parker.com

Parker's Motion & Control Technologies

At Parker, we're guided by a relentless drive to help our customers become more productive and achieve higher levels of profitability by engineering the best systems for their requirements. It means looking at customer applications from many angles to find new ways to create value. Whatever the motion or control technology need, Parker has the experience, breadth of product and global reach to consistently delivery. No company knows more about motion and control technology than Parker. For further information call 00800 27 27 5374



AEROSPACE

Key Markets

- Aircraft engines
- Business & general aviation
- Commercial transports
- Land-based weapons systems
- Military aircraft
- Missiles & launch vehicles
- Regional transports
- Unmanned aerial vehicles

Key Products

- Flight control systems & components
- Fluid conveyance systems
- Fluid metering delivery & atomization devices
- Fuel systems & components
- Hydraulic systems & components
- Inert nitrogen generating systems
- Pneumatic systems & components
- Wheels & brakes



CLIMATE CONTROL

Key Markets

- Agriculture
- Air conditioning
- Food, beverage & dairy
- Life sciences & medical
- Precision cooling
- Processing
- Transportation

Key Products

- CO₂ controls
- Electronic controllers
- Filter driers
- Hand shut-off valves
- Hose & fittings
- Pressure regulating valves
- Refrigerant distributors
- Safety relief valves
- Solenoid valves
- Thermostatic expansion valves



ELECTROMECHANICAL

Key Markets

- Aerospace
- Factory automation
- Life science & medical
- Machine tools
- Packaging machinery
- Paper machinery
- Plastics machinery & converting
- Primary metals
- Semiconductor & electronics
- Textile
- Wire & cable

Key Products

- AC / DC drives & systems
- Electric actuators, gantry robots & slides
- Electrohydraulic actuation systems
- Electromechanical actuation systems
- Human machines interface
- Linear motors
- Stepper motors, servo motors, drives & controls
- Structural extrusions



FILTRATION

Key Markets

- Food & beverage
- Industrial machinery
- Life sciences
- Marine
- Mobile equipment
- Oil & gas
- Power generation
- Process
- Transportation

Key Products

- Analytical gas generators
- Compressed air & gas filters
- Condition monitoring
- Engine, air, fuel & oil filtration & systems
- Process, chemical, water & microfiltration filters
- Nitrogen, hydrogen & zero air generators



FLUID & GAS HANDLING

Key Markets

- Aerospace
- Agriculture
- Bulk chemical handling
- Construction machinery
- Food & beverage
- Fuel & gas delivery
- Industrial machinery
- Mobile
- Oil & gas
- Transportation
- Welding

Key Products

- Brass fittings & valves
- Diagnostic equipment
- Fluid conveyance systems
- Industrial hose
- PTFE & PFA hose, tubing & plastic fittings
- Rubber & thermoplastic hose & couplings
- Tube fittings & adapters
- Quick disconnects




HYDRAULICS

Key Markets

- Aerospace
- Aerial lift
- Agriculture
- Construction machinery
- Forestry
- Industrial machinery
- Mining
- Oil & gas
- Power generation & energy
- Truck hydraulics

Key Products

- Diagnostic equipment
- Hydraulic cylinders & accumulators
- Hydraulic motors & pumps
- Hydraulic systems
- Hydraulic valves & controls
- Power take-offs
- Rubber & thermoplastic hose & couplings
- Tube fittings & adapters
- Quick disconnects



PNEUMATICS

Key Markets

- Aerospace
- Conveyor & material handling
- Factory automation
- Life science & medical
- Machine tools
- Packaging machinery
- Transportation & automotive

Key Products

- Air preparation
- Brass fittings & valves
- Manifolds
- Pneumatic accessories
- Pneumatic actuators & grippers
- Pneumatic valves & controls
- Quick disconnects
- Rotary actuators
- Rubber & thermoplastic hose & couplings
- Structural extrusions
- Thermoplastic tubing & fittings
- Vacuum generators, cups & sensors



PROCESS CONTROL

Key Markets

- Chemical & refining
- Food, beverage & dairy
- Medical & dental
- Microelectronics
- Oil & gas
- Power generation

Key Products

- Analytical sample conditioning products & systems
- Fluoropolymer chemical delivery fittings, valves & pumps
- High purity gas delivery fittings, valves & regulators
- Instrumentation fittings, valves & regulators
- Medium pressure fittings & valves
- Process control manifolds



SEALING & SHIELDING

Key Markets

- Aerospace
- Chemical processing
- Consumer
- Energy, oil & gas
- Fluid power
- General industrial
- Information technology
- Life sciences
- Military
- Semiconductor
- Telecommunications
- Transportation

Key Products

- Analytical sample conditioning products & systems
- Dynamic seals
- Elastomeric o-rings
- EMI shielding
- Extruded & precision-cut, fabricated elastomeric seals
- Homogeneous & inserted elastomeric shapes
- High temperature metal seals
- Metal & plastic retained composite seals
- Thermal management



ENGINEERING YOUR SUCCESS.

Parker Worldwide

AE – UAE, Dubai
Tel: +971 4 8127100
parker.me@parker.com

AR – Argentina, Buenos Aires
Tel: +54 3327 44 4129

AT – Austria, Wiener Neustadt
Tel: +43 (0)2622 23501-0
parker.austria@parker.com

AT – Eastern Europe, Wiener Neustadt
Tel: +43 (0)2622 23501 900
parker.easteurope@parker.com

AU – Australia, Castle Hill
Tel: +61 (0)2-9634 7777

AZ – Azerbaijan, Baku
Tel: +994 50 2233 458
parker.azerbaijan@parker.com

BE/LU – Belgium, Nivelles
Tel: +32 (0)67 280 900
parker.belgium@parker.com

BR – Brazil, Cachoeirinha RS
Tel: +55 51 3470 9144

BY – Belarus, Minsk
Tel: +375 17 209 9399
parker.belarus@parker.com

CA – Canada, Milton, Ontario
Tel: +1 905 693 3000

CH – Switzerland, Etoy
Tel: +41 (0) 21 821 02 30
parker.switzerland@parker.com

CL – Chile, Santiago
Tel: +56 2 623 1216

CN – China, Shanghai
Tel: +86 21 2899 5000

CZ – Czech Republic, Klecany
Tel: +420 284 083 111
parker.czechrepublic@parker.com

DE – Germany, Kaarst
Tel: +49 (0)2131 4016 0
parker.germany@parker.com

DK – Denmark, Ballerup
Tel: +45 43 56 04 00
parker.denmark@parker.com

ES – Spain, Madrid
Tel: +34 902 330 001
parker.spain@parker.com

FI – Finland, Vantaa
Tel: +358 (0)20 753 2500
parker.finland@parker.com

FR – France, Contamine s/Arve
Tel: +33 (0)4 50 25 80 25
parker.france@parker.com

GR – Greece, Athens
Tel: +30 210 933 6450
parker.greece@parker.com

HK – Hong Kong
Tel: +852 2428 8008

HU – Hungary, Budapest
Tel: +36 1 220 4155
parker.hungary@parker.com

IE – Ireland, Dublin
Tel: +353 (0)1 466 6370
parker.ireland@parker.com

IN – India, Mumbai
Tel: +91 22 6513 7081-85

IT – Italy, Corsico (MI)
Tel: +39 02 45 19 21
parker.italy@parker.com

JP – Japan, Tokyo
Tel: +[81] 3 6408 3901

KR – South Korea, Seoul
Tel: +82 2 559 0400

KZ – Kazakhstan, Almaty
Tel: +7 7272 505 800
parker.easteurope@parker.com

LV – Latvia, Riga
Tel: +371 6 745 2601
parker.latvia@parker.com

MX – Mexico, Apodaca
Tel: +52 81 8156 6000

MY – Malaysia, Shah Alam
Tel: +60 3 7849 0800

NL – The Netherlands, Oldenzaal
Tel: +31 (0)541 585 000
parker.nl@parker.com

NO – Norway, Ski
Tel: +47 64 91 10 00
parker.norway@parker.com

NZ – New Zealand, Mt Wellington
Tel: +64 9 574 1744

PL – Poland, Warsaw
Tel: +48 (0)22 573 24 00
parker.poland@parker.com

PT – Portugal, Leca da Palmeira
Tel: +351 22 999 7360
parker.portugal@parker.com

RO – Romania, Bucharest
Tel: +40 21 252 1382
parker.romania@parker.com

RU – Russia, Moscow
Tel: +7 495 645-2156
parker.russia@parker.com

SE – Sweden, Spånga
Tel: +46 (0)8 59 79 50 00
parker.sweden@parker.com

SG – Singapore
Tel: +65 6887 6300

SK – Slovakia, Banská Bystrica
Tel: +421 484 162 252
parker.slovakia@parker.com

SL – Slovenia, Novo Mesto
Tel: +386 7 337 6650
parker.slovenia@parker.com

TH – Thailand, Bangkok
Tel: +662 717 8140

TR – Turkey, Istanbul
Tel: +90 216 4997081
parker.turkey@parker.com

TW – Taiwan, Taipei
Tel: +886 2 2298 8987

UA – Ukraine, Kiev
Tel: +380 44 494 2731
parker.ukraine@parker.com

UK – United Kingdom, Warwick
Tel: +44 (0)1926 317 878
parker.uk@parker.com

US – USA, Cleveland
Tel: +1 216 896 3000

VE – Venezuela, Caracas
Tel: +58 212 238 5422

ZA – South Africa, Kempton Park
Tel: +27 (0)11 961 0700
parker.southafrica@parker.com

European Product Information Centre
Free phone: 00 800 27 27 5374
(from AT, BE, CH, CZ, DE, DK, EE, ES, FI, FR, IE, IL, IS, IT, LU, MT, NL, NO, PL, PT, RU, SE, UK, ZA)

© 2009 Parker Hannifin Corporation

07/09 Rev. 1A

Parker Hannifin Ltd
Parker domnick hunter
Process Division
Durham Road
Birtley, Co. Durham
DH3 2SF, England
phone +44 (0)191 4105121
fax +44 (0)191 4105312
email: dhprocess@parker.com
www.domnickhunter.com

Parker Process Advanced Filtration
2340 Eastman Avenue
Oxnard, CA 93030, USA
phone +1 805 604 3400
fax +1 805 604 3401

